


SPOT SAMPLER™ Aerosol Particle Collector



Collects concentrated samples of nanometer to micrometer-sized aerosol particles for chemical or biological analysis.

Technology: Water condensation enlarges particle size for gentle inertial deposition into a solid or liquid collection terminus; high-quality samples are ready for analysis with minimal handling.

- Applications:***
- Environmental Air Quality and Health Effects
 - Bioaerosol Collection
 - Semiconductor Micro-particle Contamination



*Collection substrates:
a) multi-well plate; b) liquid vial; and c) silicon wafer on SEM stub.*



Developed by Aerosol Dynamics Inc. with NIEHS SBIR support

Licensed and manufactured by Aerosol Devices Inc., Fort Collins, CO
For more information: www.aerosoldevices.com
www.aerosol.us